

RaphCare mobile update (v1.8.48+60)

Android patient app build **1.8.48** (install code **60**). File name: `RaphCare-v1.8.48+60.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Makes the vendor-scan crash clue readable. Build 1.8.47 could wipe the step file when scan stopped (before Connect), and the Loading overlay only said "Loading..." with no step code.

Fixed

- Stopping the watch SDK scan no longer clears the crash-step file before Connect.
- After a crash, Devices or Watch readings should show the step again (including SCAN-STOP).
- Scan setup errors are caught and recorded as SCAN-PROXY-FAILED instead of only "keeps stopping".

Added

- Live step text on the Loading overlay during **Try vendor scan (diagnostic)** (for example Vendor probe: SCAN-INVOKE).

Updated

- Version label **1.8.48** and install code **60**.
- Same single-stack vendor scan probe. Keep **1.8.42** for normal Connect.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.48+60.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → **Try vendor scan (diagnostic)**.
2. Watch the Loading text for step codes (SCAN-1, SCAN-INVOKE, SCAN-RESULT, SCAN-STOP, CONNECT-...).
3. If the app closes, reopen Devices or Watch readings and read the red step message.
4. Tell your engineer the last live step and the relaunch step (they should match or be the next step).
5. Optional: try the official H Band / Veepoo sample app on this phone with a clean reboot between tests.