

RaphCare mobile update (v1.8.16+28)

Android patient app build **1.8.16** (install code **28**). File name: `RaphCare-v1.8.16+28.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Tapping Measure now on Watch readings could close the whole app.

Fixed

- Measure no longer waits for the watch connect reply on the phone UI thread (that deadlock looked like the app crashing).
- Measure still retries when the watch radio is busy, with a clear message if it cannot connect.

Updated

- Version label **1.8.16** and install code **28**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.16+28.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Connect your claimed band.
2. Open Watch readings. On the watch open Heart Rate and tap to test.
3. Tap Measure now. The app should stay open. You should get a heart rate, a clear busy-radio message, or a stop message, not a close to the home screen.
4. Switch to Home. Content should still be there.