

RaphCare mobile update (v1.8.19+31)

Android patient app build **1.8.19** (install code **31**). File name: `RaphCare-v1.8.19+31.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Measure was still force-closing the app while the watch was measuring. The watch SDK path is unsafe on these sample bands right now, so this build stops calling it.

Fixed

- Measure no longer uses the watch vendor SDK path that was closing RaphCare (with or without Tap to test on the watch).
- The app should stay open when you tap Measure. You will see a clear on-screen note that live Measure is temporarily off.
- Connect on Devices is unchanged. The Bluetooth link should still stay up when you leave the page.

Updated

- Version label **1.8.19** and install code **31**.
- Watch readings explains that Measure is paused while we fix the crash.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.19+31.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Connect your claimed band. Confirm it stays Connected when you leave and return.
2. Open Watch readings and tap Measure now. The app must stay open and show that live Measure is temporarily off (or a reading if the phone Bluetooth link already has one).
3. Do not expect heart rate from Measure on this build. That is intentional until the watch SDK crash is fixed.