

RaphCare mobile update (v1.8.20+32)

Android patient app build **1.8.20** (install code **32**). File name: `RaphCare-v1.8.20+32.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Measure kept closing the app because the phone used two Bluetooth stacks at once. This build connects the watch through the watch SDK only, then Measure starts heart rate on that same session.

Fixed

- Connect no longer leaves a normal Bluetooth link that Measure later steals for the watch SDK (that handoff was crashing the app).
- Measure only runs after a successful watch SDK Connect, and only asks for heart rate (not oxygen at the same time).
- The app stays on the watch SDK session after Measure instead of switching Bluetooth stacks again.

Updated

- Version label **1.8.20** and install code **32**.
- Watch readings reminds you to stay on the screen while measuring.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.20+32.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Scan, then Connect your claimed band. Wait until Connected succeeds.
2. Open Watch readings. Keep RaphCare in the foreground.
3. On the watch open Heart Rate and tap to test if needed.
4. Tap Measure now. You should get a heart rate, or a clear on-screen message. The app must not close by itself.
5. If it still closes, tell us. With USB debugging on, we can capture a phone log with `scripts/Capture-RaphCareAndroidLogcat.ps1`.