

RaphCare mobile update (v1.8.22+34)

Android patient app build **1.8.22** (install code **34**). File name: `RaphCare-v1.8.22+34.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: 1.8.21 could sit on "Reconnecting your claimed watch..." and show a red Scan message without ever finishing. This build either reconnects, or it stops and lets you Scan.

Fixed

- Opening Devices no longer stays on Reconnecting. If the watch is nearby it connects. If not, status tells you to Scan, then Connect.
- Reconnect no longer depends only on extra watch libraries on the phone. If those are missing, Devices does a short Scan and Connects the claimed band itself.

Updated

- Version label **1.8.22** and install code **34**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.22+34.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build over 1.8.21.
2. Open Devices. You may see Reconnecting for about 10 seconds. It should then show Connected, or a clear line to Scan then Connect. It must not stay on Reconnecting.
3. If it asks you to Scan, tap Scan, then Connect on the claimed band (ET585).
4. Close RaphCare fully, open Devices again, and confirm the same: Connected, or Scan then Connect. Not a stuck Reconnecting line.
5. Open Watch readings and Measure. The app must not close by itself.