



RaphCare mobile update (v1.8.23+35)

Android patient app build **1.8.23** (install code **35**). File name: `RaphCare-v1.8.23+35.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Scan failed with a raw Nearby devices permission error after the app tried to reconnect the claimed watch.

Fixed

- Scan, Connect, and auto-reconnect ask for Nearby devices permission on the phone UI thread before talking to Bluetooth.
- If permission is still missing, you see a clear Settings message instead of a long system error about BLUETOOTH_SCAN.

Updated

- Version label **1.8.23** and install code **35**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.23+35.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build. When Android asks for Nearby devices (or Bluetooth), tap Allow.
2. Open Devices and tap Scan. You should not see the long BLUETOOTH_SCAN system error.
3. If permission was denied earlier, open phone Settings, find RaphCare, turn on Nearby devices, return to the app, and Scan again.
4. Connect the claimed band, then try Watch readings Measure.