

RaphCare mobile update (v1.8.24+36)

Android patient app build **1.8.24** (install code **36**). File name: `RaphCare-v1.8.24+36.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: After auto-connect showed Connected, Measure still said the watch was not ready for readings, and left "Measuring..." on screen with the red error.

Fixed

- Measure reuses the watch link from auto-connect, or opens it once if the phone only had a normal Bluetooth Connected state.
- When Measure fails, the grey status line no longer stays stuck on "Measuring..."

Updated

- Version label **1.8.24** and install code **36**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.24+36.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install and allow Nearby devices if Android asks.
2. Open Devices. Let the claimed band auto-connect, or Scan and Connect.
3. Open Watch readings. Keep RaphCare in the foreground.
4. On the watch open Heart Rate and tap to test if needed.
5. Tap Measure now. You should get a heart rate, or a clear message without a stuck "Measuring..." line. The app must not close by itself.