

RaphCare mobile update (v1.8.25+37)

Android patient app build **1.8.25** (install code **37**). File name: `RaphCare-v1.8.25+37.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Measure was closing the app (sometimes with "RaphCare keeps stopping"). The watch SDK was told to stop a heart reading that had never been started.

Fixed

- Measure no longer calls stop-reading on the watch SDK before the first start-reading.
- Measure no longer tries a second Bluetooth handshake mid-Measure (that handoff was also closing the app).

Updated

- Version label **1.8.25** and install code **37**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.25+37.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Force-stop the separate **H Band** app if it is installed (phone Settings, Apps), so it is not holding the watch.
2. Install this build. Allow Nearby devices if Android asks.
3. Devices: wait for Connected (auto-connect or Scan then Connect).
4. Watch readings → Measure now. You can open Heart Rate on the watch first, or skip that if a reading is already on the watch face.
5. The app must stay open. You should get a heart rate, or a clear on-screen message, not a system crash dialog.