

RaphCare mobile update (v1.8.27+39)

Android patient app build **1.8.27** (install code **39**). File name: `RaphCare-v1.8.27+39.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Auto-connect worked, but tapping Scan showed an empty Nearby list. The phone was still holding the watch link, so Scan could not see it, and the name filter could hide the claimed band.

Fixed

- Manual Scan releases an existing watch link first so the band can advertise again.
- Filtered Scan always includes your claimed Bluetooth address, even when the watch name does not look like E580 / E585.

Updated

- Version label **1.8.27** and install code **39**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.27+39.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build.
2. Open Devices. If it auto-connects, that is fine.
3. Tap Scan. Nearby should list your claimed band (status may briefly leave Connected while Scan runs).
4. Tap Connect, then Watch readings → Measure now.