

RaphCare mobile update (v1.8.47+59)

Android patient app build **1.8.47** (install code **59**). File name: `RaphCare-v1.8.47+59.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Fixes the engineer **Try vendor scan (diagnostic)** button on Devices. On 1.8.46 the button showed but stayed disabled after the claimed watch loaded.

Fixed

- **Try vendor scan (diagnostic)** enables after Devices loads your claimed watch (and refreshes correctly when busy state changes).

Updated

- Version label **1.8.47** and install code **59**.
- Same vendor-scan probe as 1.8.46. Keep **1.8.42** for normal Connect.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.47+59.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Open Devices. Confirm **Try vendor scan (diagnostic)** looks tappable (not greyed out).
2. Tap it. Status should say the vendor scan probe is running.
3. If the app closes, reopen Devices or Watch readings and read the red step code.
4. If Connect survives, open Watch readings and try Measure.