

RaphCare mobile update (v1.8.17+29)

Android patient app build **1.8.17** (install code **29**). File name: `RaphCare-v1.8.17+29.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Measure now could show "This feature is not supported" and then force-close RaphCare.

Fixed

- Measure no longer calls watch disconnect APIs before the watch SDK has an active session (that path caused the unsupported toast and the crash).
- Measure first tries the existing Bluetooth link for a few seconds before switching to the watch SDK.
- Safer method matching when talking to the watch SDK, so the wrong Bluetooth call is less likely.

Updated

- Version label **1.8.17** and install code **29**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.17+29.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Connect your claimed band.
2. Open Watch readings. On the watch open Heart Rate and tap to test.
3. Tap Measure now. The app must stay open. You should get a reading, or a clear on-screen message, not "RaphCare keeps stopping".
4. You should not see a lasting "This feature is not supported" toast that closes the app.