

RaphCare mobile update (v1.8.18+30)

Android patient app build **1.8.18** (install code **30**). File name: `RaphCare-v1.8.18+30.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Measure could run for a while (watch measuring), then close RaphCare with no message.

Fixed

- Measure now takes heart rate only. It no longer starts blood oxygen while heart measure is still running (that overlap closed the app).
- Measure stops the watch detect cleanly before reconnecting normal Bluetooth.
- Safer handling of watch SDK callbacks so a bad reply cannot crash the process as easily.

Updated

- Version label **1.8.18** and install code **30**.
- Blood oxygen on Measure stays for a later build; heart rate is the focus here.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.18+30.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Connect your claimed band.
2. Open Watch readings. On the watch open Heart Rate and tap to test.
3. Tap Measure now. Keep the app in the foreground. You should get a heart rate, or a clear on-screen message. The app must not close by itself.
4. Blood oxygen may stay blank on this build. That is expected for now.