

RaphCare mobile update (v1.8.21+33)

Android patient app build **1.8.21** (install code **33**). File name: `RaphCare-v1.8.21+33.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: closing the app dropped the watch Bluetooth session. Opening Devices again and tapping Connect could close the app. This build reconnects your claimed watch when you come back.

Fixed

- Closing RaphCare no longer leaves you stuck disconnected. Open Devices again and the claimed watch reconnects when it is nearby. You should not need to Scan first.
- Tapping Connect on a watch that is already connected no longer tears the session down and rebuilds it (that could close the app).
- Watch readings Measure still works after that reconnect.

Updated

- Version label **1.8.21** and install code **33**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.21+33.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Devices: Connect your claimed band (ET585 in the last test). Wait until the row says Connected.
2. Close RaphCare fully (swipe it away from recent apps).
3. Open RaphCare, then Devices. Status should return to Connected while the band is nearby. You should not have to Scan first.
4. The other nearby band (for example ET580) should still be rejected if it is not the claimed watch.
5. Open Watch readings. On the watch open Heart Rate and tap to test if needed. Tap Measure now. You should get a heart rate, or a clear message. The app must not close by itself.