

RaphCare mobile update (v1.8.30+42)

Android patient app build **1.8.30** (install code **42**). File name: `RaphCare-v1.8.30+42.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: On 1.8.29, Scan could finish with an empty Nearby list and no Claimed wearable row. That happens when the claimed device has no Bluetooth address from the server. This build reloads the claim before Scan and shows a clear message instead of a blank list.

Fixed

- Scan reloads the claimed device and checks for a usable Bluetooth MAC before searching.
- If the MAC is missing, you see an explicit message (not a silent empty Nearby list).
- Debug logs record the API MAC and normalized MAC for support.

Updated

- Version label **1.8.30** and install code **42**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.30+42.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build.
2. Devices → Scan.
3. If you see a red message about a missing Bluetooth address: that is the real issue. Ask Ops to set the MAC on the device, or use Show all Bluetooth, Connect once to lock it, then Scan again.
4. If the MAC is present, you should see **Claimed wearable** (or the band), then Connect and Measure.