

RaphCare mobile update (v1.8.42+54)

Android patient app build **1.8.42** (install code **54**). File name: `RaphCare-v1.8.42+54.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Probe 1.8.41 turned watch-SDK Connect back on and Scan → Connect still closed the app. This build returns to the proven phone Bluetooth Connect path.

Fixed

- Connect uses phone Bluetooth again so the app stays open after Scan.
- Watch-SDK Connect and live Measure stay off until a safe native path exists.

Updated

- Version label **1.8.42** and install code **54**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.42+54.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Watch nearby, unlocked, on your wrist.
2. Devices → Scan until your ET585 appears.
3. Tap Connect. The app must stay open and reach Connected.
4. Do not expect Measure through the watch SDK in this build.