

RaphCare mobile update (v1.8.43+55)

Android patient app build **1.8.43** (install code **55**). File name: `RaphCare-v1.8.43+55.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Isolated probe for watch-SDK Connect. Expect Scan → Connect may still close the app. Keep **1.8.42** for normal use. The point of this build is to capture a useful crash log with the fixed capture script.

Updated

- Watch-SDK Connect and live Measure are on for diagnosis only.
- Version label **1.8.43** and install code **55**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.43+55.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.
5. Keep `RaphCare-v1.8.42+54.apk` ready to reinstall after the test.

What to try (with log capture)

1. On the PC, with the phone on USB debugging: `powershell -NoProfile -ExecutionPolicy Bypass -File scripts/Capture-RaphCareAndroidLogcat.ps1`
2. Start the script **before** you tap Connect.
3. Watch nearby, unlocked, on your wrist.
4. Devices → Scan until ET585 appears → tap Connect.
5. If the app stays open, try Watch readings → Measure.
6. Share the file under `artifacts/android/logs/` (look for CONNECT-2 without CONNECT-3, Fatal signal, DEBUG tombstone).
7. If Connect closed the app, reinstall **1.8.42**.