

RaphCare mobile update (v1.8.29+41)

Android patient app build **1.8.29** (install code **41**). File name: `RaphCare-v1.8.29+41.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Manual Scan could finish with an empty Nearby list even when the claimed watch was already known. Auto-connect by address worked; the scan list did not.

Fixed

- Nearby keeps a **Claimed wearable** row for your locked Bluetooth address when Scan finds no advertisement.
- Connect on that row uses the same watch SDK path as auto-reconnect.
- The row shows as soon as Scan starts, not only after the long scan window.
- Devices no longer stays on Connected after a failed watch SDK handshake.

Updated

- Version label **1.8.29** and install code **41**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.29+41.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build. Force-stop the separate **H Band** app if installed.
2. Devices → Scan. You should see **Claimed wearable** even if no other rows appear.
3. Tap Connect on that row. Wait until Connected succeeds.
4. Watch readings → Measure now.