

RaphCare mobile update (v1.8.41+53)

Android patient app build **1.8.41** (install code **53**). File name: `RaphCare-v1.8.41+53.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Probe for watch-SDK Connect and live Measure after init fixes. Connect may still close the app on some phones. Keep **1.8.40** nearby if you need the stable Connect path.

Updated

- Watch-SDK Connect and live Measure are on again for this probe.
- Opening Devices warms up the watch SDK without connecting yet.
- After Scan, the app waits briefly before watch-SDK Connect so the radio can settle.
- Version label **1.8.41** and install code **53**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.41+53.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Watch nearby, unlocked, on your wrist.
2. Open Devices (wait a moment for warm-up).
3. Scan until your ET585 appears, then tap Connect.
4. The app must stay open and reach Connected.
5. Open Watch readings and try Measure. Look for a heart rate sample.
6. If Connect closes the app, reinstall `RaphCare-v1.8.40+52.apk` and tell us.