

RaphCare mobile update (v1.8.49+61)

Android patient app build **1.8.49** (install code **61**). File name: `RaphCare-v1.8.49+61.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: USB log capture is unreliable on this Huawei. The app now writes a probe log on the phone you can Share (email, Drive, WhatsApp) without adb.

Added

- **Share probe log** on Devices (sends `raphcare-vendor-probe-log.txt`).
- **Test breadcrumb (no watch)** to confirm the step file works before a real scan.
- Every vendor step is appended to that log (it is not wiped when Connect clears the single-step breadcrumb).

Updated

- Version label **1.8.49** and install code **61**.
- Same vendor-scan diagnostic as 1.8.48.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.49+61.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → **Test breadcrumb (no watch)**. You should see Breadcrumb OK and a file path.
2. Tap **Share probe log** and send the file to yourself (or save to Drive / Files).
3. Then **Try vendor scan (diagnostic)**. After a crash, reopen the app and tap **Share probe log** again (even if there is no red step).
4. Also note the last Loading text (`Vendor probe: ...`) if you saw it.

Where the file lives on the phone

`Android/data/com.yindula.raphcare/files/raphcare-vendor-probe-log.txt`

(Share is easier than digging that folder on Huawei.)