

RaphCare mobile update (v1.8.50+62)

Android patient app build **1.8.50** (install code **62**). File name: `RaphCare-v1.8.50+62.apk` . This build talks to the staging API, not a phone on localhost.

Why this build: The on-phone probe log from 1.8.49 showed vendor scan and connect both finished, then the app died before handshake. This build records the next steps (wait for callbacks, password, person sync) and keeps a red step after relaunch if connect returned but handshake did not.

Fixed

- After a crash past connect, Devices can show the last incomplete step again (connect return no longer clears that marker).

Added

- Finer probe steps after connect: wait for Bluetooth connect, wait for notify, password, person sync.
- Same **Share probe log** and **Test breadcrumb (no watch)** as 1.8.49.

Updated

- Version label **1.8.50** and install code **62**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.50+62.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → **Try vendor scan (diagnostic)**.
2. If the app closes, reopen it. Note any red step text (for example WAIT-CONNECT or PWD-1).
3. Tap **Share probe log** and send `raphcare-vendor-probe-log.txt` .
4. Also note the last Loading text if you saw it.

Where the file lives on the phone

`Android/data/com.yindula.raphcare/files/raphcare-vendor-probe-log.txt`

(Share is easier than digging that folder on Huawei.)