

RaphCare mobile update (v1.8.51+63)

Android patient app build **1.8.51** (install code **63**). File name: `RaphCare-v1.8.51+63.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: The probe log from 1.8.50 stopped at WAIT-CONNECT (connect started, then the app closed while waiting for the watch Bluetooth callback). This build waits briefly after stopping the scan before connecting, and holds the connect listeners more firmly.

Fixed

- Short pause after vendor scan stop before connect (same settle used elsewhere).
- Connect and notify listeners stay pinned until handshake finishes.

Updated

- Version label **1.8.51** and install code **63**.
- Same Share probe log and Test breadcrumb tools as before.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.51+63.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → **Try vendor scan (diagnostic)**. Loading may pause about 1.5 seconds after the watch is found.
2. If it still closes, reopen and note the red step (SETTLE-1, WAIT-CONNECT, WAIT-NOTIFY, PWD-1, or HANDSHAKE if you get further).
3. Tap **Share probe log** and send the file.

Where the file lives on the phone

`Android/data/com.yindula.raphcare/files/raphcare-vendor-probe-log.txt`