

RaphCare mobile update (v1.8.44+56)

Android patient app build **1.8.44** (install code **56**). File name: `RaphCare-v1.8.44+56.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Probe for watch-SDK Connect when USB logging is not available. If Connect closes the app, reopen Devices and read the red message (it names the last step, for example CONNECT-2). Keep **1.8.42** for normal use.

Added

- After a Connect crash, Devices can show the last watch-SDK step so you can report it without a PC log.

Updated

- Watch-SDK Connect and Measure stay on for this probe only.
- Version label **1.8.44** and install code **56**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.44+56.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.
5. Keep `RaphCare-v1.8.42+54.apk` ready to reinstall after the test.

What to try

1. Watch nearby, unlocked, on your wrist.
2. Devices → Scan → Connect.
3. If the app closes, open RaphCare again and open Devices. Read the red error (note the step code like CONNECT-2).
4. Tell your engineer that step code, then reinstall **1.8.42**.
5. Optional: if USB debugging works, also run the capture script before Connect.