

RaphCare mobile update (v1.8.45+57)

Android patient app build **1.8.45** (install code **57**). File name: `RaphCare-v1.8.45+57.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Fixes the crash-step message. Build 1.8.44 wrote the step, then auto-reconnect cleared it and often crashed again before you could read it.

Fixed

- After a Connect crash, Devices shows the step code and pauses auto-reconnect so the red message stays on screen.

Updated

- Version label **1.8.45** and install code **57**.
- Still a watch-SDK Connect probe. Keep **1.8.42** for normal use.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.45+57.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → Scan → Connect.
2. If the app closes, open RaphCare and open Devices once.
3. Read the red message (look for CONNECT-2 or CONNECT-INVOKE). Do not tap Connect again until you have noted the code.
4. Tell your engineer that code, then reinstall **1.8.42**.