

RaphCare mobile update (v1.9.5+71)

Android patient app build **1.9.5** (install code **71**). File name: `RaphCare-v1.9.5+71.apk` . This build talks to the staging API, not a phone on localhost.

Why this build: On Huawei, vendor scan Connect stayed open through password start (1.9.4), then failed with “confirmDevicePwd overload not found.” The watch SDK method uses a primitive boolean; our call used a boxed Boolean and reflection never matched. This build fixes that match.

Keep for daily Connect: `RaphCare-v1.8.42+54.apk` .

Fixed

- Password confirm after vendor Connect now finds the watch SDK method (boolean argument match).

Updated

- Version label **1.9.5** and install code **71**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Install `RaphCare-v1.9.5+71.apk` .
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.
5. Prefer forgetting the ET585 in phone Bluetooth settings before the first vendor scan.

What to try

1. Devices → **Try vendor scan (diagnostic)** with the watch nearby and unlocked.
2. Ideal: trail reaches PERSON-1 / HANDSHAKE-OK without “confirmDevicePwd overload not found.”
3. Then Watch readings → Measure.
4. If anything fails, tap **Share probe log** and send `raphcare-vendor-probe-log.txt` .