

RaphCare mobile update (v1.9.6+72)

Android patient app build **1.9.6** (install code **72**). File name: `RaphCare-v1.9.6+72.apk` . Staging API.

Why this build: Vendor Connect handshake now works (HANDSHAKE-OK on Huawei). Measure was taking an early “still detecting” heart sample, so the phone could show a lower number than the watch face. Watch readings also showed “Heart rate” twice.

Keep for daily Connect: `RaphCare-v1.8.42+54.apk` .

Fixed

- Measure waits for a settled heart reading (not the interim detecting sample).
- Watch readings shows the bpm number once under the Heart rate label.

Updated

- Version label **1.9.6** and install code **72**.

Install

1. Install `RaphCare-v1.9.6+72.apk` (uninstall older RaphCare only if Android blocks the update).
2. Force-stop the separate **H Band** app if it is installed.
3. Sign in as the demo patient.

What to try

1. Devices → **Try vendor scan (diagnostic)** until Connect succeeds.
2. Watch readings → wear the band snug → **Measure now** (do not rely on the watch’s own Heart Rate screen alone).
3. Expect a bpm close to what the watch shows after it finishes measuring.
4. If wrong, **Share probe log** from Devices.