

RaphCare mobile update (v1.8.52+64)

Android patient app build **1.8.52** (install code **64**). File name: `RaphCare-v1.8.52+64.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: 1.8.51 still stopped at WAIT-CONNECT even after a pause after scan stop. This build keeps the watch SDK scan running while Connect starts, then stops the scan afterward.

Updated

- Vendor diagnostic: find watch, connect while scan is still on, stop scan in cleanup.
- Probe log shows **SCAN-KEEP** before Connect steps.
- Version label **1.8.52** and install code **64**.

Before you try (phone settings)

1. Android Bluetooth settings: forget the ET585 if it is listed, then toggle Bluetooth off and on.
2. Huawei battery / App launch: set RaphCare to manual or unrestricted (not auto-manage).

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.52+64.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Devices → **Try vendor scan (diagnostic)**.
2. If it closes, reopen and note the red step. Share the probe log.
3. Ideal progress: past WAIT-CONNECT to WAIT-NOTIFY or a successful Connect.