

RaphCare mobile update (v1.8.40+52)

Android patient app build **1.8.40** (install code **52**). File name: `RaphCare-v1.8.40+52.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Turning watch-SDK Connect back on (even with the documented address-only call) still closed the app after Scan. This build returns to the proven phone Bluetooth Connect path from 1.8.38.

Fixed

- Connect uses phone Bluetooth again so the app stays open after Scan.
- Watch-SDK Connect and live Measure stay off until a safe native path exists.

Updated

- Version label **1.8.40** and install code **52**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.40+52.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Watch nearby, unlocked, on your wrist.
2. Devices → Scan until your ET585 appears.
3. Tap Connect. The app must stay open and reach Connected.
4. Do not expect Measure through the watch SDK in this build.