

RaphCare mobile update (v1.8.28+40)

Android patient app build **1.8.28** (install code **40**). File name: `RaphCare-v1.8.28+40.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Scan still showed no Nearby watches. The phone was dropping blank-named bands before it could match your claimed Bluetooth address, and paired watches were not seeded into the list.

Fixed

- Manual Scan keeps advertisements long enough to match your claimed address (including blank watch names).
- Already paired or system-connected watches that match your claim appear in Nearby.
- Clearer empty-scan hint (force-stop H Band, try Show all Bluetooth).

Updated

- Version label **1.8.28** and install code **40**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.28+40.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Install this build. Force-stop the separate **H Band** app if it is on the phone.
2. Devices → Scan. Nearby should list your claimed band (or a paired match).
3. If still empty: tap the Showing row to list all Bluetooth, Scan again, then Connect your claimed band.
4. Watch readings → Measure now.