

RaphCare mobile update (v1.8.46+58)

Android patient app build **1.8.46** (install code **58**). File name: `RaphCare-v1.8.46+58.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Two engineer fixes for watch Connect diagnosis. Normal Connect stays on the stable **1.8.42** path. The crash-step message now shows on Watch readings as well as Devices. A separate diagnostic button runs watch-SDK scan then connect without the usual Bluetooth scan stack.

Fixed

- After a Connect crash, Watch readings also shows the step code and pauses auto-reconnect (same as Devices).
- Crash-step handling is shared so a later screen cannot skip it by mistake.

Added

- Devices button **Try vendor scan (diagnostic)** for engineer tests only. Uses the watch SDK for scan and connect in one path.

Updated

- Version label **1.8.46** and install code **58**.
- Daily Connect remains Plugin.BLE style (**1.8.42** behavior). Keep **1.8.42** installed if you need normal Connect while testing.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.46+58.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Claim the watch if needed. Do **not** use the normal Scan button for this probe.
2. Devices → **Try vendor scan (diagnostic)**.
3. If the app closes, open RaphCare (Devices or Watch readings). Read the red message (SCAN-... or CONNECT-...). Note the code before tapping Connect again.
4. If Connect survives, open Watch readings and try Measure.
5. Tell your engineer the result, then reinstall **1.8.42** for daily use.