

RaphCare mobile update (v1.9.0+66)

Android patient app build **1.9.0** (install code **66**). File name: `RaphCare-v1.9.0+66.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Preview checkpoint for the commercial wave (Practice plan, consult waiting display, pricing, and payments work that follows). This is **not** the daily stable Connect APK. Keep using your last good Connect build for watch tests until we say otherwise.

Added

- Version label **1.9.0** and install code **66** as the commercial-wave preview marker.

Updated

- Same Connect path as 1.8.53 for now (mac-only overload and warm scan). No new watch fix in this preview alone.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.9.0+66.apk` to the phone.
3. Force-stop the separate **H Band** app if it is installed.
4. Sign in as the demo patient.

What to try

1. Sign in and open Home to confirm the build label.
2. For watch Connect diagnostics, prefer your last known-good Connect APK unless you are specifically testing this preview.