

RaphCare mobile update (v1.8.26+38)

Android patient app build **1.8.26** (install code **38**). File name: `RaphCare-v1.8.26+38.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: Devices showed Connected while Watch readings Measure still said to connect first. The phone had a normal Bluetooth link, not the measure link Measure needs.

Fixed

- Opening Devices or Watch readings refreshes the measure link when Connected is only Bluetooth.
- Connect no longer leaves a Bluetooth-only Connected state when the watch SDK is available.
- Clearer on-screen guidance when Measure still needs Disconnect then Connect.

Updated

- Version label **1.8.26** and install code **38**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.26+38.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Force-stop the separate **H Band** app if installed.
2. Install this build.
3. Open Devices. Wait for Connected (or Scan then Connect). Status should stay Connected after the measure link is ready.
4. Open Watch readings → Measure now.
5. If Measure still asks you to reconnect: on Devices tap Disconnect, then Connect, then Measure again. The app must stay open.