

RaphCare partner update (6 Sep 2026)

Here is what landed today. The short version: patients can **measure heart rate** from a claimed E580 or E585 band on Watch readings, Ops **stocks a watch from the packaging first** (serial and Bluetooth address on the confirm step), and the chat assistant **keeps the last few turns** in the same session. Current Android build is **1.8.15**. Overall partner checklist completion is about **83%**.

1. Patient Android app (1.8.15)

Sideload file name: `RaphCare-v1.8.15+27.apk`. It talks to the staging API, not localhost. Copy it to the phone (USB, Drive, or email). Allow install from that source if Android asks. If an older RaphCare install blocks it, uninstall the old copy first. This build is for sideload testing, not the Play Store.

What changed since the 1.8.4 note:

- **Watch readings.** After claim and Connect, open Watch readings and tap Measure now. Heart rate (and oxygen when the band sends it) can fill on the phone. Using Health Monitor on the watch alone does not send numbers. Measure retries if the watch radio is still busy, and stops with a clear message instead of spinning forever.
- **Connection that lasts.** Leaving Devices keeps the Bluetooth link when the band is still nearby. Connect only accepts the claimed watch after the first pair. A wrong nearby band is rejected with the address that belongs to your claim.
- **Claim from a photo.** On Devices, photograph the packaging barcode or QR, then Claim watch. You can still type the serial.
- **Emergency contacts from the phone book.** On Profile, pick a contact or type one in.
- **Chat that remembers this session.** The assistant can use the last few turns, so a follow-up makes sense. This is still not a doctor.
- **Calmer look, fewer crashes.** Main screens look more open. Notifications, Devices, Active clinic, Health records scroll, and opening AI from Home should stay up instead of closing the app.

One-page phone note: [mobile-update-v1.8.15+27.pdf](#).

2. Ops fleet (packaging first)

On the Ops site, adding a watch is now two steps. Capture the packaging first (barcode, picture in the browser, or a photo you already have). Then confirm hospital, serial, Bluetooth address, and model. The picture stays in the browser. It is not uploaded.

Bluetooth address is required. The fleet list shows it when known. Ops can still assign, take a watch back to stock, or remove a mistaken row (history is retired, not erased). Hospital staff still assign an in-stock watch to a patient on Portal Devices.

3. Feature catalog

The shareable list of what RaphCare does is updated for Watch readings, claim-and-lock pairing, and the Ops capture-first stock flow: [raphcare-feature-catalog-v1.pdf](#). It is not a progress board. Completion % stays on the partner status sheet.

4. Demo clinic and accounts (unchanged)

Staging still seeds **RaphCare Demo Clinic** when the API starts. It does not wipe other hospitals.

Shared password for every demo account: **RaphCareDemo!2026**

Who	Email
Hospital admin	demo.admin@raphcare.com
Doctor	demo.doctor@raphcare.com
Pharmacist	demo.pharmacy@raphcare.com
Lab	demo.lab@raphcare.com
Platform Ops	demo.ops@raphcare.com
Patient (phone app)	demo.patient@raphcare.com

Staff sign-in: [Professional sign-in](#)

Ops sign-in: [Ops sign-in](#)

A one-page sheet with the same logins is here: [raphcare-demo-accounts-v2.pdf](#). Keep that sheet private. It has the password.

5. Links

- Portal: <https://raphcare-hqf6gsa3acanargz.southafricanorth-01.azurewebsites.net>
- Ops: <https://raphcare-ops.azurewebsites.net>
- API: <https://raphcare-api-eydjcnefhae2dpa2.southafricanorth-01.azurewebsites.net>
- Public site: raphcare.com

6. Partner checklist

The partner status board is updated for this work. Rough rollup:

Area	%
Hospital admin (outpatient)	100%
Inpatient	100%
Patient phone app	61%
Other staff and system	93%
Hospital plan extras (stay and counter)	100%
Overall	83%

Try this week

1. Install `RaphCare-v1.8.15+27.apk` and sign in as `demo.patient@raphcare.com`.
2. Open Devices. Claim the assigned watch if needed (photo of the packaging code, or type the serial), then Scan and Connect. Leave Devices and come back: status should still say Connected while the band is nearby.
3. Open Watch readings. On the watch, open Heart Rate and tap to test. Tap Measure now. You should get a heart rate, or a clear retry message (not an endless spinner).
4. Open the chat assistant. Send a short wellness question, then a follow-up that only makes sense with the first message. The reply should still know what you were talking about.
5. Open Profile, Emergency contacts. Pick someone from the phone book (or type a name).
6. On Ops, sign in as `demo.ops@raphcare.com`. Add a watch: capture packaging first, then confirm serial and Bluetooth address.