

RaphCare mobile update (v1.8.7+19)

Android patient app build **1.8.7** (install code **19**). File name: `RaphCare-v1.8.7+19.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: claim an assigned watch by photographing the packaging barcode or QR, not only by typing the serial.

Added

- On Devices and wearables, **Photo of barcode** fills the claim serial from a clear packaging barcode or QR (then tap Claim watch).

Updated

- Claim hint explains photo or typed serial.
- Version label **1.8.7** and install code **19**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.7+19.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Sign in as the demo patient.
2. Open Devices and wearables.
3. Tap **Photo of barcode**, photograph the packaging code for an assigned watch, confirm the serial field fills, then Claim watch.
4. You can still type the serial if the photo does not read cleanly.