



RaphCare mobile update (v1.8.8+20)

Android patient app build **1.8.8** (install code **20**). File name: `RaphCare-v1.8.8+20.apk`. This build talks to the staging API, not a phone on localhost.

Why this build: clearer layout on Devices, Request a call, and Mental health so buttons stay on screen and the screens feel calmer.

Updated

- **Devices and wearables:** card layout, clearer Scan / Claim / Connect actions, softer inputs.
- **Request a call:** Video, Audio, and Low data sit in one equal-width control so nothing is cut off.
- **Mental health:** mood choices in a two-by-two grid; PHQ-9 and GAD-7 stack full width.
- Version label **1.8.8** and install code **20**.

Install

1. Uninstall any older RaphCare build if Android blocks the update.
2. Copy `RaphCare-v1.8.8+20.apk` to the phone (USB, Drive, or email).
3. Allow install from that source if Android asks.
4. Sign in as the demo patient: `demo.patient@raphcare.com` (password on the demo accounts sheet).

This is for sideload testing, not the Play Store.

What to try

1. Sign in as the demo patient.
2. Open Devices and wearables. Check that Scan, Stop scan, and Claim fit the screen.
3. Open Request a call. Confirm Video, Audio, and Low data are all visible.
4. Open Mental health. Confirm mood buttons and Start PHQ-9 / Start GAD-7 are fully on screen.